



 **SolM** Building
a Healthier
Tomorrow

USA CATALOG

VETERINARY CARE PRODUCTS



solm.com | info@solm.com

©2023 Sol-Millennium Medical Inc. All rights reserved.

Sol-Millennium is a vertically integrated global medical supply company with locations across the globe. We are one of the world's largest manufacturers of needles and syringes and have a comprehensive line of medical supplies.



Our innovative technology reduces needlestick injuries, saves on medication expenses and enhances both clinician and patient experience. Our international sales team is dedicated to delivering knowledge and outstanding customer service. Our goal is to work together as partners in the healthcare industry, focusing on a healthier tomorrow.



- Global HQ
- Research & Development
- Sales, Marketing & Distribution
- Manufacturing

▶ Table of Contents



	Needles	3-6
	Syringes	7-10
	Diabetes Care	11-12
	IV Infusion	13-14
	Blood Collection	15-18
	Protective Apparel	19-22
	Gloves	23-24

Sol - Vet™ products are for veterinary use only.

Sol - Vet™

● New Vet Grind

Our commitment to excellence give you world class quality, reliability, precision, and performance in every box. Available in most popular gauges and lengths

Features:

- Triple-bevel lubricated stainless steel bevel that minimizes patient discomfort
- Point durability optimized for vet procedure
- Color-coded hubs for easy needle gauge identification
- Available in soft packs



Vet Grind


Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
V111810	Sol-Vet Hypodermic Needle 18G x 1"	18G	1"	100	1,000
V111815	Sol-Vet Hypodermic Needle 18G x 1-1/2"	18G	1-1/2"	100	1,000
V112110	Sol-Vet Hypodermic Needle 20G x 1"	20G	1"	100	1,000
V112015	Sol-Vet Hypodermic Needle 20G x 1-1/2"	20G	1-1/2"	100	1,000
V112034	Sol-Vet Hypodermic Needle 20G x 3/4"	20G	3/4"	100	1,000
V112210	Sol-Vet Hypodermic Needle 22G x 1"	22G	1"	100	1,000
V112234	Sol-Vet Hypodermic Needle 22G x 3/4"	22G	3/4"	100	1,000
V112506	Sol-Vet Hypodermic Needle 25G x 5/8"	25G	5/8"	100	1,000

● Hypodermic Needles

Sol-Millennium® offers a wide range of hypodermic needles with an emphasis on ease of use and patient comfort.

Features:

- Triple-bevel lubricated stainless steel bevel that minimizes patient discomfort
- Color-coded hubs for easy gauge identification
- Compatible with all standard Luer Lock and Slip Tip syringes
- Available in a wide range of gauges and needle lengths



Hypodermic Needles

Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
111910	Sol-M Hypodermic Needle 19G x 1"	19G	1"	100	1,000
111915	Sol-M Hypodermic Needle 19G x 1-1/2"	19G	1-1/2"	100	1,000
112110	Sol-M Hypodermic Needle 21G x 1"	21G	1"	100	1,000
112115	Sol-M Hypodermic Needle 21G x 1-1/2"	21G	1-1/2"	100	1,000
112215	Sol-M Hypodermic Needle 22G x 1 1/2"	22G	1-1/2"	100	1,000
112310	Sol-M Hypodermic Needle 23G x 1"	23G	1"	100	1,000
112358	Sol-M Hypodermic Needle 23G x 5/8"	23G	5/8"	100	1,000
112506	Sol-M Hypodermic Needle 25G x 5/8"	25G	5/8"	100	1,000
112515	Sol-M Hypodermic Needle 25G x 1-1/2"	25G	1-1/2"	100	1,000
112705	Sol-M Hypodermic Needle 27G x 1/2"	27G	1/2"	100	1,000
113005	Sol-M Hypodermic Needle 30G x 1/2"	30G	1/2"	100	1,000



NEEDLES


Sol-Care™

● Safety Needles

Sol-Care™ Safety Needles help to prevent needlestick injuries. Every needle is equipped safety shield mechanism that locks in place to lower the risk of accidental injury after medication administration.

Features:

- Audible click indicates safety mechanism is secure and locked in place
- Protective arm does not obstruct visualization of injection site
- Bevel-up needle allows for low-angle injections
- Compatible with all standard Luer Lock and Luer Slip Tip syringes



Sol-Care™ Safety Needles

Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
SN1810	Sol-Care Safety Needle 18G x 1"	18G	1"	100	800
SN1815	Sol-Care Safety Needle 18G x 1-1/2"	18G	1-1/2"	100	800
SN2010	Sol-Care Safety Needle 20G x 1"	20G	1"	100	800
SN2015	Sol-Care Safety Needle 20G x 1-1/2"	20G	1-1/2"	100	800
SN2110	Sol-Care Safety Needle 21G x 1"	21G	1"	100	800
SN2210	Sol-Care Safety Needle 22G x 1"	22G	1"	100	800
SN2215	Sol-Care Safety Needle 22G x 1-1/2"	22G	1-1/2"	100	800
SN2310	Sol-Care Safety Needle 23G x 1"	23G	1"	100	800
SN2510	Sol-Care Safety Needle 25G x 1"	25G	1"	100	800
SN2558	Sol-Care Safety Needle 25G x 5/8"	25G	5/8"	100	800
SN2515	Sol-Care Safety Needle 25G x 1-1/2"	25G	1-1/2"	100	800
SN2712	Sol-Care Safety Needle 27G x 1/2"	27G	1/2"	100	800
SN2758	Sol-Care Safety Needle 27G x 5/8"	27G	5/8"	100	800
SN3012	Sol-Care Safety Needle 30G x 1/2"	30G	1/2"	100	800

Sol-Vet™

Needle and Syringe Combos

Sol-Vet Needle and Syringe Combos are designed to match the specific needs of the Veterinarian.

Features:

- Combines our Veterinary engineered needles with high quality syringes
- Syringes have wide finger flanges for easier control and use
- Available in most popular luer slip and luer lock configurations
- Available in soft packs for reduced waste



Sol-Vet™ Needle and Syringe Combos

Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
VSSC001	Sol-Vet 1mL Luer Slip w/detachable 25G x 5/8"	25G	5/8"	100	800
VSSC009	Sol-Vet 3mL Luer Slip w/detachable 20G x 1"	20G	1"	100	800
VSSC010	Sol-Vet 3mL Luer Slip w/detachable 20G x 1 1/2"	20G	1 1/2"	100	800
VSSC011	Sol-Vet 3mL Luer Slip w/detachable 22G x 3/4"	22G	3/4"	100	800
VSSC012	Sol-Vet 3mL Luer Slip w/detachable 22G x 1"	22G	1"	100	800
VSSC013	Sol-Vet 3mL Luer Slip w/detachable 22G x 1 1/2"	22G	1 1/2"	100	800
VSSC002	Sol-Vet 3mL Luer Lock w/20G x 1"	20G	1"	100	800
VSSC003	Sol-Vet 3mL Luer Lock w/20G x 1 1/2"	20G	1 1/2"	100	800
VSSC004	Sol-Vet 3mL Luer Lock w/20G x 3/4"	20G	3/4"	100	800
VSSC005	Sol-Vet 3mL Luer Lock w/22G x 3/4"	22G	3/4"	100	800
VSSC006	Sol-Vet 3mL Luer Lock w/22G x 1"	22G	1"	100	800
VSSC007	Sol-Vet 3mL Luer Lock w/23G x 1"	23G	1"	100	800
VSSC013	Sol-Vet 3mL Luer Lock w/25G x 5/8"	25G	5/8"	100	800



NEEDLES


Sol-M™

3-Piece Syringes

Sol-M™ 3-Piece Syringes accurately measure the quantity of liquid medications and are available in a range of sizes and connection combinations to address a variety of clinical needs.

Features:

- High-transparency barrel to easily visualize contents
- Silicone lubrication ensures smooth plunger movement
- Plunger backstop that prevents accidental withdrawal of plunger from the barrel
- Large finger grips allow for easy, controlled injections
- Available in Luer Lock and Slip Tip options



Luer Lock Syringe
without needle

Luer Slip Tip Syringe
without needle



Sol-M™ 3-Piece Syringes

Ref. Number	Product Description	Volume	Box Qty	Case Qty
*P180001PP	Sol-M 1mL Luer Lock Syringe w/o Needle (PP) featuring Low Dead Space (LDS) technology	1 mL	100	800
*P180001LDS	Sol-M 1mL Luer Lock Syringe w/o Needle featuring Low Dead Space (LDS) technology	1mL	100	800
P180003	Sol-M 3mL Luer Lock Syringe w/o Needle	3 mL	100	800
180003ST	Sol-M 3mL Slip Tip Syringe w/o Needle	3 mL	100	800
P180005	Sol-M 5mL Luer Lock Syringe w/o Needle	5 mL	100	800
180005ST	Sol-M 5mL Slip Tip Syringe w/o Needle	5 mL	100	800
P180010	Sol-M 10mL Luer Lock Syringe w/o Needle	10 mL	100	1,200
180010ST	Sol-M 10mL Slip Tip Syringe w/o Needle	10 mL	100	1,200
P180020	Sol-M 20mL Luer Lock Syringe w/o Needle	20 mL	100	800
180020ST	Sol-M 20mL Slip Tip Syringe w/o Needle	20 mL	100	800
P180030	Sol-M 30mL Luer Lock Syringe w/o Needle	30 mL	50	400
P180060	Sol-M 60mL Luer Lock Syringe w/o Needle	60 mL	30	240

*Low Dead Space (LDS) syringes are designed to minimize the amount of residual medication waste left within the syringe after medication has been delivered.

Sol-M™

Oral Dispensing Syringes

The Sol-M™ Oral Syringe has been specifically designed for use with oral medications, utilizing a special tip that will not accept Luer connectors, which lowers the risk of medication errors.

Features:

- 3-piece syringe with O-ring and gasket
- Large variety of sizes and colors
- Clear graduation markings (mL, tsp) for accurate medication delivery
- Optional amber barrel designed for light-sensitive drugs, in compliance with USP requirements
- Protective tip cap for transportation with no spillage



Available with clear barrel or with an amber barrel designed for light-sensitive drugs

Sol-M™ Oral Dispensing Syringes – Gasket Type - Bulk Non-Sterile

Ref. Number	Product Description	Barrel Color	Volume	Box Qty	Case Qty
21010GCNSB	Sol-M 10mL Oral Dispensing Syringe Clear w/Tip Cap (Gasket, Bulk-Non-Sterile)	Clear	1 mL	100	400
21010GANSB	Sol-M 1mL Oral Dispensing Syringe Amber w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Amber	1 mL	100	400
21005GCNSB	Sol-M 5mL Oral Dispensing Syringe Clear w/Tip Cap (Gasket, Bulk-Non-Sterile)	Clear	5 mL	100	400
21005GANSB	Sol-M 5mL Oral Dispensing Syringe Amber w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Amber	5 mL	100	400
21003GCNSB	Sol-M 3mL Oral Dispensing Syringe Clear w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Clear	3 mL	100	400
21003GANSB	Sol-M 3mL Oral Dispensing Syringe Amber w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Amber	3 mL	100	400
21001GANSB	Sol-M 1mL Oral Dispensing Syringe Amber w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Clear	1 mL	100	400
21001GCNSB	Sol-M 1mL Oral Dispensing Syringe Clear w/ Tip Cap (Gasket, Bulk-Non-Sterile)	Amber	1 mL	100	400



Sol-M™ Oral Dispensing Syringes – O-ring Type

Ref. Number	Product Description	Barrel Color	Volume	Box Qty	Case Qty
21001OA	Sol-M 1mL Oral Dispensing Syringe Amber w/ Tip Cap	Amber	1 mL	100	400
21001OC	Sol-M 1mL Oral Dispensing Syringe Clear w/ Tip Cap	Clear	1 mL	100	400
21003OA	Sol-M 3mL Oral Dispensing Syringe Amber w/ Tip Cap	Amber	3 mL	100	400
21003OC	Sol-M 3mL Oral Dispensing Syringe Clear w/ Tip Cap	Clear	3 mL	100	400
21005OA	Sol-M 5mL Oral Dispensing Syringe Amber w/ Tip Cap	Amber	5 mL	100	400
21005OC	Sol-M 5mL Oral Dispensing Syringe Clear w/ Tip Cap	Clear	5 mL	100	400
20010OA	Sol-M 10mL Oral Dispensing Syringe Amber w/ Tip Cap	Amber	10 mL	100	400
21010OC	Sol-M 10mL Oral Dispensing Syringe Clear w/ Tip Cap	Clear	10 mL	100	400

Sol-M™ Oral Dispensing Syringes – Gasket Type

Ref. Number	Product Description	Barrel Color	Volume	Box Qty	Case Qty
21001GA	Sol-M 1mL Oral Dispensing Syringe Amber w/Tip Cap (Gasket type)	Amber	1 mL	100	800
21001GC	Sol-M 1mL Oral Dispensing Syringe Clear w/Tip Cap (Gasket type)	Clear	1 mL	100	800
21003GA	Sol-M 3mL Oral Dispensing Syringe Amber w/Tip Cap (Gasket type)	Amber	3 mL	100	800
21003GC	Sol-M 3mL Oral Dispensing Syringe Clear w/Tip Cap (Gasket type)	Clear	3 mL	100	800




● Insulin Syringes

Sol-Millennium's® Insulin Syringes provide accurate and comfortable insulin administration.

Features:

- Needles engineered to ensure comfortable insulin administration
- Bold and clear graduated insulin unit markings for accurate dosing
- Available in multiple sizes



(U-40)



Insulin Syringes U-40

Ref. Number	Product Description	Volume	Needle Length	Qty. per bag	Qty. bags per box	Qty. bags per case
V3305-40	Sol-Vet 0.3ml U-40 Insulin Syringe w/Fixed Needle 30G x 1/2"(half unit markings)(PE Bag)	0.3 mL	1/2"	10	10	100
V32905-40	Sol-Vet 0.3ml U-40 Insulin Syringe w/Fixed Needle 29G x 1/2"(half unit markings)(PE Bag)	0.3 mL	1/2"	10	10	100
V5295-40	Sol-Vet 0.5mL U-40 Insulin Syringe w/Fixed Needle 29G x 1/2"(half unit markings)(PE Bag)	0.5 mL	1/2"	10	10	100
V5305-40	Sol-Vet 0.5mL U-40 Insulin Syringe w/Fixed Needle 30G x 1/2"(half unit markings)(PE Bag)	0.5 mL	1/2"	10	10	100
V12905-40	Sol-Vet 1mL U-40 Insulin Syringe w/Fixed Needle 29G x 1/2" (PE Bag)	1 mL	1/2"	10	10	100
V1305-40	Sol-Vet 1mL U-100 Insulin Syringe w/Fixed Needle 30G x 1/2" (PE Bag)	1 mL	1/2"	10	10	100

Insulin Syringes U-100

Ref. Number	Product Description	Volume	Needle Length	Qty. per bag	Qty. bags per box	Qty. bags per case
V3295	Sol-Vet 0.3mL U-100 Insulin Syringe w/Fixed Needle 29G x 1/2" (half unit markings)(PE Bag)	0.3 mL	1/2"	10	10	100
V3305	Sol-Vet 0.3mL U-100 Insulin Syringe w/Fixed Needle 30G x 1/2" (half unit markings)(PE Bag)	0.3 mL	1/2"	10	10	100
V5295	Sol-Vet 0.5mL U-100 Insulin Syringe w/Fixed Needle 29G x 1/2" (half unit markings)(PE Bag)	0.5 mL	1/2"	10	10	100
V5305	Sol-Vet 0.5mL U-100 Insulin Syringe w/Fixed Needle 30G x 1/2" (half unit markings) (PE Bag)	0.5 mL	1/2"	10	10	100
V12905	Sol-Vet 1mL U-100 Insulin Syringe w/Fixed Needle 29G x 1/2" (PE Bag)	1 mL	1/2"	10	10	100
V1305	Sol-Vet 1mL U-100 Insulin Syringe w/Fixed Needle 30G x 1/2" (half unit markings)(PE Bag)	1mL	1/2"	10	10	100


Sol-Care™

● Safety Winged Infusion Sets

Sol-Care™ Safety Winged Infusion Sets are intended for short-term IV medication administration. Each infusion set is equipped with a safety mechanism that is activated after the infusion has completed, helping to prevent needlestick injuries.

Features:

- Thin-walled needle minimizes patient discomfort
- DEHP-free tubing is soft, flexible and kink-resistant



Sol-Care™ Safety Winged Infusion Set

Ref. Number	Product Description	Gauge	Needle Length	Tubing Length	Case Qty
213412SI	Sol-Care Safety Winged Infusion Set 21G x 3/4 x 12"	21G	3/4"	12"	200
253412SI	Sol-Care Safety Winged Infusion Set 25G x 3/4 x 12"	25G	3/4"	12"	200

Sol-Vet™

● IV Catheters

Sol-Vet™ Safety IV Catheters were designed for veterinarians to have a safety feature on the IV medication administration for animals.

Features:

- Protective casing ensures safety of catheter during transportation
- Translucent catheter hub and chamber allows for easy visualization of blood flashback
- Contains a hydrophobic membrane filter to prevent blood leakage during insertion
- Filter cap can be removed to allow connection with a syringe to assist insertion



Sol-Vet™ IV Catheters

Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
110301050001	Sol-Vet IV Catheter 18G x 1-1/4"	18G	1-1/4"	100	1,000
110301050002	Sol-Vet IV Catheter 18G x 1-3/4"	18G	1-3/4"	100	1,000
110301050007	Sol-Vet IV Catheter 20G x 1"	20G	1"	100	1,000
110301050003	Sol-Vet IV Catheter 20G x 1-1/4"	20G	1-1/4"	100	1,000
110301050004	Sol-Vet IV Catheter 22G x 1"	22G	1"	100	1,000
110301050005	Sol-Vet IV Catheter 24G x 3/4"	24G	3/4"	100	1,000
110301050006	Sol-Vet IV Catheter 26G x 3/4"	26G	3/4"	100	1,000





Sol - Guard™



● Safety Pull-Button Blood Collection Sets

Sol-Guard™ Safety Pull-Button Blood Collection Sets promote single-handed in-vein safety activation to instantly retract the needle after use, helping prevent needlestick injury and bloodborne pathogen exposure.

With ease, safety and comfort in mind, this device features an exclusive version of the “Butterfly Needle”, designed with a sliding button that using pull force, activates in-vein needle retraction after use.


Features:

- Spring mechanism for split second in-vein needle retraction
- Large oversized arrow to easily identify the activation pull button
- Flashback for instant and clear confirmation of proper needle placement
- Low angled tri-bevel tip designed to improve penetration, while minimizing patient discomfort
- Thin-wall needle for fast flow rate and to reduce the risk of blood hemolysis
- Tactile non-slip grip to help engage the safety feature
- 30% wider wings for improved handling during insertion and removal



Available with or without pre-attached holder

Sol-Guard™ Safety Pull-Button Blood Collection Set with Pre-attached Holder



Blood flashback and flow are easy to visualize

Sol-Guard™ Safety Pull-Button Blood Collection Set

Providing a **Safe, Easy & Comfortable** Blood Collection Experience



BLOOD COLLECTION

88%

Reduction in needlestick injuries with in-vein needle retraction safety activation¹


2,000,000

Estimated number of needlestick and other sharps related injuries annually²

3

41%

of phlebotomist incur injuries from blood-draw procedures annually³



Sol-Guard® Safety Pull-Button Blood Collection Sets

Reference number		Needle outer diameter		Needle length		Tubing length	
without pre-attached holder	with pre-attached holder	Gauge	Millimeters	Inches	Millimeter	Inches	Millimeter
-	PBLA2107	21G	0.8 mm	3/4"	19 mm	7"	178 mm
-	PBLA2307	23G	0.6 mm	3/4"	19 mm	7"	178 mm
-	PBLA2507	25G	0.5 mm	3/4"	19 mm	7"	178 mm
PBH2112	PBLA2112	21G	0.8 mm	3/4"	19 mm	12"	305 mm
PBH2312	PBLA2312	23G	0.6 mm	3/4"	19 mm	12"	305 mm
PBH2512	PBLA2512	25G	0.5 mm	3/4"	19 mm	12"	305 mm

Sales Unit With pre-attached holder: 25 pcs/Box; 200 pcs/Case
Without pre-attached holder: 50 pcs/Box; 200 pcs/Case



Sol-M™

Multi-Sample Needles

The Sol-M™ Multi-Sample Needle makes it possible to collect multiple samples with a single vein puncture.

Features:

- Marking on cap to indicate which direction provides bevel up orientation
- Low-angle bevel design allows easy penetration
- Color-coded caps and hubs allow easy identification of needle gauge size
- Compatible with all standard holders



Black dot indicates bevel up orientation of needle

Color-coded caps and hubs allow easy identification of needle gauge size



Low-angle bevel design allows easy penetration



Sol-M™ Multi-Sample Needles

Ref. Number	Product Description	Gauge	Needle Length	Box Qty	Case Qty
110201020008	Sol-M Multi-Sample Needle 18G x 1"	18G	1"	100	1,000
110201020009	Sol-M Multi-Sample Needle 18G x 1-1/2"	18G	1-1/2"	100	1,000




● Butterfly Needles

The Sol-Vet™ Butterfly Needle features a removeable multi-sample luer adapter, making its application flexible across multiple uses.

The Sol-Care™ Safety Blood Collection Needle is a safety-engineered device with a sliding safety shield that clicks when activated, telling you the needle is safely locked.

Features:

- Thin-walled needle minimizes patient discomfort
- Transparent design allows for clear visualization of blood flow
- Available as a safety engineered device (Sol-Care™ Safety Blood Collection Needle)
- Available with removeable multi-sample luer adapter (Sol-Vet™ Butterfly Needle)
- Available with or without pre-attached blood collection tube holder



Butterfly Needles

Ref. Number	Product Description	Gauge	Needle Length	Tubing Length	Case Qty
110301070004	Sol-Vet Butterfly Needle with Adapter 19G x 3/4 x 12"	19G	3/4"	12"	1,000
110301070001	Sol-Vet Butterfly Needle with Adapter 21G x 3/4 x 12"	21G	3/4"	12"	1,000
110301070002	Sol-Vet Butterfly Needle with Adapter 23G x 3/4 x 12"	23G	3/4"	12"	1,000
110301070003	Sol-Vet Butterfly Needle with Adapter 25G x 3/4 x 12"	25G	3/4"	12"	1,000

Blood Collection Accessories

Ref. Number	Product Description	Box Qty	Case Qty
110201050001	Sol-M Blood Collection Tube Holder	250	1,000
110201050003	Sol-M Holder with Male Luer Lock Adapter	200	800



Beard Cover

Ref. Number	Description	Box	Case
SGCWB01	Sol-M Beard Cover	100	1,000

Surgeons Cap with Ties

Ref. Number	Description	Box	Case
SGCBT01	Sol-M Surgeons Cap w/ Ties	100	600

Protective Face Shield

Ref. Number	Description	Box	Case
FS9132	Sol-M Protective Face Shield	1	300

Surgeons Hood

Ref. Number	Description	Box	Case
SGCWH01	Sol-M Surgeons Hood	100	500

Bouffant Caps

Ref. Number	Description	Size	Box	Case
BFCB21	Sol-M Bouffant Cap – Blue	21"	100	1,000
BFCB24	Sol-M Bouffant Cap – Blue	24"	100	1,000
BFCB27	Sol-M Bouffant Cap – Blue	27"	100	1,000
BFCW21	Sol-M Bouffant Cap – White	21"	100	1,000
BFCW24	Sol-M Bouffant Cap – White	24"	100	1,000
BFCW27	Sol-M Bouffant Cap – White	27"	100	1,000

Face Masks

Ref. Number	Product Description	Box Qty	Case Qty
MSK048	Sol-M Dinosaur Children's Masks w/ Earloops, Fold Mask	10	500
KMSK001	Sol-M Children's Protective Mask, Yellow	10	600
KMSK002	Sol-M Children's Protective Mask, Green	10	600
KMSK00150	Sol-M ASTM Level 2 Children's Mask, Earloops, Yellow	50	2,000
KMSK00250	Sol-M ASTM Level 2 Children's Mask, Earloops, Green	50	2,000
MSK003X	Sol-M ASTM Level 1 Medical Mask w/ Earloops	50	2,000
MSK003X2	Sol-M ASTM Level 2 Surgical Mask w/ Earloops	50	2,000
MSK032	Sol-M ASTM Level 2 Surgical Mask w/ Ties	50	1,000
MSK034	Sol-M ASTM Level 2 Surgical Mask w/ Eye Shield, Earloops	50	500
MSK036	Sol-M ASTM Level 2 Surgical Mask w/ Eye Shield, w/Ties	50	500
MSK033	Sol-M ASTM Level 3 Surgical Mask w/ Eye Shield, Earloops	50	500
MSK035	Sol-M ASTM Level 3 Surgical Mask w/ Eye Shield, w/Ties	50	500
MSKEA30W	Sol-M ASTM Level 3 Surgical Mask w/ Earloops	50	2,000
MSKET30W	Sol-M ASTM Level 3 Surgical Mask w/ Ties	50	1,000
MSK003KU	Sol-M Protective Mask w/ Earloops	50	2,000
MSK003XR	Sol-M Protective Face Mask w/Earloops	20	2,000

Respirators

Ref. Number	Description	Box	Case
YQD8008	Ychita Niosh N95 Respirator Cup-Shape	15	270
YQ950101	Sol- M KN95 Respirator w/ Earloops	25	900



PROTECTIVE APPAREL

Gowns

Ref. Number	Product Description	Size	Case Qty
IG0049	Sol-M Chemotherapy Gown, Blue	Large	80
IG0050	Sol-M Chemotherapy Gown, Blue	XLarge	50
IG232U1111	Sol-M AAMI Level 1 Isolation Gown, Ties, Yellow	Universal	100
IG232XL1111	Sol-M AAMI Level 1 Isolation Gown, Ties, Yellow	XLarge	100
IG004	Sol-M AAMI Level 2 Isolation Gown Over-the-Head, White	Large	100
IG004XL	Sol-M AAMI Level 2 Isolation Gown Over-the-Head, White	XLarge	100
IG004XXL	Sol-M AAMI Level 2 Isolation Gown Over-the-Head, White	XXLarge	100
IG0035	Sol-M AAMI Level 3 Reinforced Surgical Gown/Towel, Blue	Medium	50
IG433U2133	Sol-M AAMI Level 3 Isolation Gown, Blue	Universal	50
IG433XL2133	Sol-M AAMI Level 3 Isolation Gown, Blue	XLarge	50
SIGA321	Sol-M AAMI Level 3 Surgical Gown, Blue	Large	50
SIGA322	Sol-M AAMI Level 3 Surgical Gown, Blue	XLarge	50
SIGA323	Sol-M AAMI Level 3 Surgical Gown, Blue	XXLarge	50
IG0036	Sol-M AAMI Level 4 Reinforced Surgical Gown/Towel, Blue	Medium	50
SIGA411	Sol-M AAMI Level 4 Reinforced Surgical Gown, Blue	Large	50
SIGA412	Sol-M AAMI Level 4 Reinforced Surgical Gown, Blue	XLarge	50
SIGA413	Sol-M AAMI Level 4 Reinforced Surgical Gown, Blue	XXLarge	50



Lab Jackets

Ref. Number	Description	Size	Bag	Case
LAC001S	Sol-M 3-Pocket Lab Coat, White	Small	10	100
LAC001M	Sol-M 3-Pocket Lab Coat, White	Medium	10	100
LAC001L	Sol-M 3-Pocket Lab Coat, White	Large	10	100
LAC001XL	Sol-M 3-Pocket Lab Coat, White	XLarge	10	100
LAC001XXL	Sol-M 3-Pocket Lab Coat, White	XXLarge	10	100
LAC0013XL	Sol-M 3-Pocket Lab Coat, White	3XLarge	10	100
LAC0014XL	Sol-M 3-Pocket Lab Coat, White	4XLarge	10	100
LAC002S	Sol-M 3-Pocket Lab Coat, Blue	Small	10	100
LAC002M	Sol-M 3-Pocket Lab Coat, Blue	Medium	10	100
LAC002L	Sol-M 3-Pocket Lab Coat, Blue	Large	10	100
LAC002XL	Sol-M 3-Pocket Lab Coat, Blue	XLarge	10	100
LAC002XXL	Sol-M 3-Pocket Lab Coat, Blue	XXLarge	10	100
LAC0023XL	Sol-M 3-Pocket Lab Coat, Blue	3XLarge	10	100
LAC0024XL	Sol-M 3-Pocket Lab Coat, Blue	4XLarge	10	100



PROTECTIVE APPAREL

Shoe Covers

Ref. Number	Description	Size	Box	Case
SC002U	Sol-M Shoe Covers – Medium Weight	Universal	100	300
SC002XL	Sol-M Shoe Covers – Medium Weight	XLarge	100	200
SC002XXL	Sol-M Shoe Covers – Medium Weight	XXLarge	100	200

Boot Covers

Ref. Number	Description	Size	Box	Case
SC004U	Sol-M Boot Covers	Universal	50	150
SC004XL	Sol-M Boot Covers	X-Large	50	150



Tested for use with Chemotherapy Drugs according to ASTM D6978-5

Sol-Millennium's nitrile examination gloves are durable enough to withstand exposure to numerous chemotherapy drugs while maintaining excellent touch-sensitivity and dexterity for clinicians.

Please see the table below for a list of all the chemotherapy drugs that have been tested with Sol-Millennium's nitrile examination gloves.

These chemotherapy drugs and concentrations did **not** have breakthrough detected up to minutes:

Tested Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Minutes)
Cisplatin (1.0 mg/mL)	No breakthrough up to 240 minutes
Cyclophosphamide (20 mg/mL)	No breakthrough up to 240 minutes
Ctarabine (100 mg/mL)	No breakthrough up to 240 minutes
Dacarbazine (10.0 mg/mL)	No breakthrough up to 240 minutes
Doxorubicin Hydrochloride (2.0 mg/mL)	No breakthrough up to 240 minutes
Etoposide (20.0 mg/mL)	No breakthrough up to 240 minutes
Fluorouracil (50.0 mg/mL)	No breakthrough up to 240 minutes
Ifosfamide (50.0 mg/mL)	No breakthrough up to 240 minutes
Methotrexate (25 mg/mL)	No breakthrough up to 240 minutes
Mitomycin C (0.5 mg/mL)	No breakthrough up to 240 minutes
Mitoxantrone (2.0 mg/mL)	No breakthrough up to 240 minutes
Paclitaxel (6.0 mg/mL)	No breakthrough up to 240 minutes
Vincristine Sulfate (1.0 mg/mL)	No breakthrough up to 240 minutes

These chemotherapy drugs and concentrations had breakthrough in less than 60 minutes

Carmustine (3.3 mg/mL)	DO NOT USE
Thiotepa (10.0 mg/mL)	DO NOT USE

WARNING: PRODUCTS NOT FOR USE WITH CARMUSTINE OR THIOTEPA

Caution: Results are based on testing conditions intended to approximate clinical usage.

It is the users' responsibility to determine if these gloves are appropriate for use with chemotherapy drugs.

● Nitrile Examination Gloves

Disposable nitrile examination gloves for medical use (ASTM D6319) that have been tested with chemotherapy drugs (ASTM D6978-5). These are powder free and not made with natural rubber latex.

Features:

- Tested for use with Chemotherapy Drugs (D6978-5)
- Powder Free
- 4.2 grams (Medium)
- Thickness: 3.5 mils
- 9.5" length
- 5-year shelf life

Ref. Number	Description	Size	Box	Case
NG55001U	SOL-M Nitrile Examination Glove	Small	100 pcs	1,000 pcs
NG55002U	SOL-M Nitrile Examination Glove	Medium	100 pcs	1,000 pcs
NG55003U	SOL-M Nitrile Examination Glove	Large	100 pcs	1,000 pcs
NG55004U	SOL-M Nitrile Examination Glove	Xlarge	100 pcs	1,000 pcs



GLOVES



solm.com | info@solm.com

